



**CUSTOM  
CONFIGURATION**

**EXPAND YOUR  
BUSINESS**

*Without Inflating Your Cost*

**Atrion**<sup>®</sup>  
Medical

Atrion Medical Products, Inc. • 1426 Curt Francis Road • Arab, AL 35016  
800 343-9334 • info@atrionmedical.com • www.atrionmedical.com

Protected under U.S. Patents 5,168,757 • 5,713,242 • 6,796,959 • 6,938,319 and DES 454,358 S. Including Foreign Patents and Patents Pending.  
QL<sup>®</sup> is a registered trademark of Atrion Medical Products, Inc.

**Atrion**<sup>®</sup>  
Medical

# Custom Built

# For Your Custom Needs

## The QL® Inflation Device

Is perfect for a wide array of applications such as balloon catheterization, stent deployment and fluid delivery.

Available in a variety of configurations so you can be sure you get exactly what you need.

## Smart Design

Its stable design and ergonomic handle make it easy to control. The quick latch mechanism allows for one-handed deflation, even at full pressure.

The clear barrel provides medical staff with a precise view of the content and volume allowing the user to observe the purging of the bubbles.

## Broad Applications

- Cardiology
- Peripheral
- Radiology
- Urological
- Gastrointestinal
- Neurology
- Orthopedics
- Biliary
- Pulmonary
- ENT

## Regulatory Approvals

Approved for use in U.S., Europe, Asia and Latin America.

## Custom Ordering

Available in a variety of configurations starting with quantities as low as 1,000. Custom name and logo can be printed on the gauge, box and label.

## Service

We work with our customers to ensure each order is built to their requirements and is delivered on time.



## Stock Configurations

| Part Number | Size | Pressure | Type        |
|-------------|------|----------|-------------|
| 96241       | 10cc | 15 ATM   | Standard    |
| 96252       | 10cc | 30 ATM   | Standard    |
| 96351       | 14cc | 55 ATM   | Low Profile |
| 96197       | 20cc | 30 ATM   | Standard    |
| 96358       | 25cc | 30 ATM   | Low Profile |
| 96311       | 25cc | 40 ATM   | Low Profile |
| 96313       | 40cc | 15 ATM   | Low Profile |
| 96314       | 60cc | 15 ATM   | Low Profile |

### QL Inflation Device At A Glance

The QL comes in two profiles:



#### Standard

Angled gauge allows for easy reading.



#### Low Profile

Easier to sterilize and easier to pack, reducing your costs.

- 510 (k) cleared
- Ergonomic handle
- Available sterile or non-sterile
- Quick latch for instant deflation
- Made in the U.S.A.
- Clear barrel
- CE marked