

Products
to advance
recovery

Cardiothoracic
& general surgery



Mission statement:

Qualiteam is dedicated to advance postoperative care and to improve the overall hospital experience for patients worldwide by developing innovative wound supports that ultimately result in cost savings for healthcare systems. We provide our products and services with a compassion as if each patient was part of our family, and we treat customers and fellow team members with honesty, integrity and respect in all aspects of our activities.

About us:

Qualiteam's production facility is located in Northern Italy, and was founded in 1999 by Doris Hjorth, a Danish cardio-thoracic nurse and a US certified perfusionist. Research and a passion to improve patient's hospital experience while enhancing the daily day patient care for nurses, are fundamental qualities rooted into the company activities.

The introduction of the patented QualiBreath sternum and thorax support in 2005 pioneered the design of external chest support devices and their use, and QualiBreath is still the only dynamic dual-functional device available.

A complete line of advanced wound support products has since been developed. The unique QualiBra Advanced sternum and breast support for fuller-sized women; the QualiBelly Advanced new generation of abdominal supports; and the postoperative care kits with exceptional products designed to help patients to a smooth recovery at home, are just examples of the growing product portfolio.

The latest addition is the Calla by Qualiteam product line, dedicated to make products for cosmetic and plastic surgery more beautiful and functionally convenient for patients and health care staff. (Not included in this catalog - see www.callabyqualiteam.com for details).

INDEX

Products - Cardiac surgery Page:

- QualiBreath sternum and thorax support - Adult
Constant lateral sternum stabilization with integrated cough support 6 - 7
- QualiBreath sternum and thorax support - Pediatric
Constant lateral sternum stabilization with integrated cough support 8 - 9
- QualiBreath sternum and thorax support - Toddler and Infant
Constant lateral sternum stabilization 10 - 11
- QualiBra Advanced sternum and breast support
Constant lateral sternum stabilization with integrated cough support and breast tissue separation 12 - 13

Products - General abdominal surgery Page:

- QualiBelly Advanced abdominal binder - Adult
Tri-band wound support with individual compression adjustment 14 - 15
- QualiBelly Advanced abdominal binder - Pediatric and Toddler
Two- and tri-band wound supports with individual compression adjustment 16 - 17
- QualiBelly post-partum, lumbar and abdominal binder
Single band wound and back support 18 - 19

| Accessories: | Page: |
|---|-------|
| • CryoPouch Non-drug adjunctive treatment for pain relief | 20 |
| • Reusable dual-purpose cool or heat packs Gel packs for cool or heat treatment | 21 |
| • AxillaPad Perspiration pads for use with QualiBreath, QualiBra Advanced or Thor products | 22 |
| • QualiPad Seatbelt pad for wound protection in the car | 23 |
| • Breast Cover For female patients with breast size below cup D wearing QualiBreath. | 24 |
| • Universal device holder and line organiser Organization of drainage systems | 25 |
| • Netting bag Washing net and container for postoperative care kit items | 26 |

| Other products: | Page: |
|--|---------|
| • QualiClip Disposable Vascular Bulldog Clip | 27 |
| • QualiMaternity belt For advanced pregnancy care | 28 - 29 |

| Gifts for patients: | Page: |
|---|---------|
| • T-shirts and scarves “Get well soon” gifts | 30 - 31 |
| • Lucky plushy toy Instructional toy dog | 32 |

QUALIBREATH

Postoperative
sternum and thorax support



- Unique nurse-designed postoperative sternum and thorax support, ideal for patients who have had a sternotomy procedure or other types of interventions/trauma in the chest area.
- Provides continuous stabilization of the sternum to prevent complications such as deep sternal wound infections and dehiscence.
- Relieves pain by a constant support to sore ribs and respiratory muscles.
- The only dual-functional sternum and thorax support:
 1. Constant, lateral sternum support which gives a uniform, encircling, elastic tension that can be adjusted.
 2. Integrated hand-grips that are always ready to reinforce the chest surrounding tension during coughing/sneezing and straining.
- Unrestricted breathing since the lower lung lobes and the upper abdominal area are free and not covered by the chest band.
- Shoulder-straps secure a correct position in bed and during physical activities.
- Available in all sizes for adults from extra small to four times extra large.
- For more information on the product features, sizes and use, please refer to the product brochure or visit www.qualiteam.com

Sizing guide and order codes:

Supplied single packed with Instructions for Use

| Measurement | Size | Order code |
|-------------------------------------|-------------|------------|
| 75 - 84 cm | Extra Small | QB XS |
| 84 - 92 cm | Small | QB S |
| 92 - 100 cm | Medium | QB M |
| 100 - 108 cm | Large | QB L |
| 108 - 113 cm | Extra Large | QB XL |
| 113 - 122 cm | XX Large | QB XXL |
| 122 - 134 cm | XXX Large | QB XXXL |
| 134 - 148 cm | 4 X Large | QB 4XL |
| Custom sizes available upon request | | |

| Accessories | |
|------------------------------------|-------------|
| Description | Order codes |
| QualiPad seatbelt pad | QTSEPA |
| AxillaPad (pair) | QBAX |
| Breast Cover (for female patients) | COV-size |

| Postoperative care kit | |
|--|--------------------|
| Postoperative care kit (QualiBreath + QualiPad + AxillaPad supplied in washing net) | QTKIT + QB size |

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.



QUALIBREATH Pediatric


Postoperative
sternum and thorax support



- Unique nurse-designed postoperative sternum and thorax support, ideal for pediatric patients who have had a sternotomy procedure or other types of interventions/trauma in the chest area.
- Provides continuous stabilization of the sternum to prevent complications and to relieve pain from sore ribs and respiratory muscles.
- The constant, lateral sternum support gives a comfortable, encircling, elastic tension that can be adjusted.
- Hand-grips integrated in the material are always in the right position for the child to grasp for reinforcing the chest surrounding tension during coughing/sneezing and when straining.
- Breathing is unrestricted since the lower lung lobes and the upper abdominal area are free and not covered by the chest band.
- Shoulder straps secure a correct position in bed and during physical activities.
- Supplied with a free plushy dog wearing a mini-version of the QualiBreath, color pencils and a coloring booklet. The booklet is composed of the "Instructions for Use" featuring the dog showing the proper use of QualiBreath to help the little patients understand how it works and protects their wound.
- For more information on the product features, sizes and use, please refer to the product brochure or visit www.qualiteam.com

Sizing guide and order codes:

Supplied single packed with Instructions for Use

|  Place measure tape ABOVE breast for correct sizing | | |
|--|----------------|------------|
| Measurement | Size | Order Code |
| 56 - 62 cm | Pediatric 4-6 | QB 4-6 |
| 62 - 68 cm | Pediatric 6-8 | QB 6-8 |
| 68 - 75 cm | Pediatric 8-10 | QB 8-10 |
| Custom sizes available upon request | | |

| Accessories | |
|--|------------|
| Description | Order Code |
| QualiPad Pediatric seatbelt pad | QTSEPAP |
| AxillaPad Pediatric (pair) | QBAXP |
| Instructional plushy dog ped. with QualiBreath | LKY-QB |

| Postoperative care kit | |
|---|--------------------|
| Postoperative care kit (QualiBreath Ped.+ QualiPad Ped.+ AxillaPad ped., Plushy dog ped. supplied in washing net) | QTKIT + QB size |



Free of charge:
Pediatric QualiBreath products are supplied with a plushy toy wearing a mini QualiBreath, along with color pencils and a coloring booklet illustrating "Instructions for Use" of QualiBreath.

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.

QUALIBREATH Infant-Toddler


Postoperative
sternum and thorax support



- Unique nurse-designed postoperative sternum and thorax support, ideal for infant and toddler patients who have had a sternotomy procedure or other types of interventions or trauma in the chest area.
- Provides continuous stabilization of the sternum to prevent complications and to relieve pain from sore ribs and respiratory muscles.
- The constant, lateral sternum support gives a comfortable, encircling, elastic tension that can be adjusted.
- Breathing is unrestricted since the lower lung lobes and the upper abdominal area are free and not covered by the chest band.
- Shoulder straps secure a correct position at all times.
- Supplied with a free plushy toy dog wearing a mini-version of the QualiBreath to help the little patients feel comforted and secure.
- For more information on the product features, sizes and use, please refer to the product brochure or visit www.qualiteam.com

Sizing guide and order codes:

Supplied single packed with Instructions for Use

|  <p>Place measure tape ABOVE breast for correct sizing</p> | | |
|---|----------------|------------|
| Measurement | Size | Order Code |
| 27 - 32 cm | Infant Small | QB INF 30 |
| 32 - 38 cm | Infant Medium | QN INF 35 |
| 38 - 46 cm | Toddler Small | QB TOD1 |
| 46 - 56 cm | Toddler Medium | QB TOD2 |
| Custom sizes available upon request | | |

| Accessories | |
|--|------------|
| Description | Order code |
| AxillaPad Infant-Toddler(pair) | QBITAX |
| Plushy dog Infant-Toddler with QualiBreath | LKY-QB |

| Postoperative care kit | |
|--|-----------------|
| Postoperative care kit (QualiBreath + AxillaPad, Plushy dog infant-toddler supplied in washing net) | QTKIT + QB size |



Infant and Toddler QualiBreath products are supplied with a plushy toy dog wearing a mini QualiBreath.

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.




Postoperative
sternum and breast support



- Unique nurse-designed postoperative sternum and breast support, ideal for the fuller-sized woman or obese patient who have had a sternotomy or other interventions in the chest.
- Provides continuous stabilization of the sternum to prevent complications such as deep sternal wound infections and dehiscence.
- Relieves pain by a constant support to sore ribs and respiratory muscles.
- Integrated pockets in the sternum band can be grasped for pain control during coughing, sneezing, straining and respiratory exercises, and provide additional support around the chest when needed.
- Wide shoulder straps and vertical breast bands can be adjusted to prevent the pull of gravity on the sternal wound, and to prevent the risk of infection by keeping it free from breast tissue to avoid heat and moisture buildup, as well as allowing for a more esthetic healing.
- Soft, elasticized cotton bra cups with no seams or labels against the skin, together with an easily adjustable front closure, provide exceptional comfort and enhance compliance for usage at home until the sternum and surgical wound has healed.
- Painless wound care while the breasts remain supported and covered.
- Available with a standard bra cup for sizes D, DD, and E, and a Plus size for cup sizes above F.
- For more information on the product features and use, please refer to the product brochure or visit www.qualiteam.com

Sizing guide and order codes:

Supplied single packed with Instructions for Use

|  Measure ABOVE breasts for correct sizing | | Standard cup size ranges D, DD, E | Plus cup size ranges F, FF and above |
|--|-------------|---|---|
| Measurement | Size | Order Code | |
| 75 - 84 cm | Extra Small | QADV-XS | QADV-XS Plus |
| 84 - 92 cm | Small | QADV-S | QADV-S Plus |
| 92 - 100 cm | Medium | QADV-M | QADV-M Plus |
| 100 - 108 cm | Large | QADV-L | QADV-L Plus |
| 108 - 113 cm | Extra Large | QADV-XL | QADV-XL Plus |
| 113 - 122 cm | XX Large | QADV-XXL | QADV-XXL Plus |
| 122 - 134 cm | XXX Large | QADV-XXXL | QADV-XXXL Plus |
| 134 - 148 cm | 4XL Large | QADV-4XL | QADV-4XL Plus |
| Custom sizes available upon request | | | |
| If the cup size is smaller than D we recommend to use the QualiBreath sternum and thorax support. | | | |

| Accessories | |
|-----------------------|------------|
| Description | Order code |
| QualiPad seatbelt pad | QTSEPA |
| AxillaPad (pair) | QBAX |

| Postoperative care kit | |
|---|----------------------|
| Postoperative care kit (QualiBra Advanced + QualiPad + AxillaPad supplied in washing net) | QTKIT + QADV size |

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.

For more information please contact Qualiteam at info@qualiteam.com
For orders email: orders@qualiteam.com

QUALIBelly advanced

Postoperative abdominal
binder



- Unique nurse-designed tri-band abdominal compression support, ideal for patients who have had major surgery in the abdominal region and have drainage systems, stoma devices, driver line systems or other devices that exit from the abdomen.
- Designed to prevent infections since separation of drainage systems, driver line systems or other devices is eliminated. The bands simply close above and below such items which exit at the slit(s) between them.
- Provides a steady pleasant support to a tender wound which reduces pain from bruising and swelling, improves blood circulation and oxygen supply to the surgical site to speed up healing of the wound.
- Individually adjusted compression settings at upper, mid and lower abdominal segments to decrease postoperative pain.
- Provides a supportive pressure to the abdominal muscles to reduce pain during increased internal abdominal pressures when coughing, sneezing and straining.
- One or two of the bands can be opened for wound access while the other band(s) remain closed to maintain support during wound care with less pain. Just one healthcare person can close the binder again on patients of all body sizes.
- Provides support to the lumbar and abdominal regions which helps to get out of bed sooner, improves posture, breathing and mobility, and to walk around more easily and safer with the hands free.
- Converts easily from a 3- to a 2-band for additional comfort when sitting.
- Convenient universal device holder available to arrange and attach drainage systems directly to the binder (sold separately).
- For more information on the product features, sizes and use, please refer to the product brochure or visit www.qualiteam.com

Sizing guide and order codes:

Supplied single packed with Instructions for Use

| | 1. Measure the fullest part of the abdomen, hips or waist. | Height: 28 cm | Height: 40 cm |
|-------------------------------------|--|--------------------------------|--------------------------------|
| | 2. Measure the distance from breast bone to pubic bone. | Measurement 2: Less than 45 cm | Measurement 2: More than 45 cm |
| Measurement 1 | Size | Order code | |
| 65 - 85 cm | Small | QLY-S | QLY40-S |
| 85 - 105 cm | Medium | QLY-M | QLY40-M |
| 105 - 125 cm | Large | QLY-L | QLY40-L |
| 125 - 145 cm | Extra Large | QLY-XL | QLY40-XL |
| 145 - 165 cm | XX Large | QLY-XXL | QLY40-XXL |
| 165 - 185 cm | XXX Large | QLY-XXXL | QLY40-XXXL |
| 185 - 205 cm | 4XL Large | QLY-4XL | QLY40-4XL |
| Custom sizes available upon request | | | |

| Accessories | |
|-------------------------|------------|
| Description | Order code |
| QualiPad seatbelt pad | QTSEPA |
| Universal device holder | DHG |

| Postoperative care kit | |
|--|---------------------------------|
| Postoperative care kit (QualiBelly Advanced + QualiPad supplied in washing net) | QTKIT + QLY or QLY40 size |

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.

For more information please contact Qualiteam at info@qualiteam.com
For orders email: orders@qualiteam.com



QUALIBelly

advanced pediatric

Postoperative abdominal
binder



- Unique nurse-designed dual- or tri-band abdominal compression support, ideal for toddler and pediatric patients who have had major surgery in the abdominal region and have drainage systems, stoma devices, driver line systems or other devices that exit from the abdomen.
- Individually adjusted compression settings at upper, mid and/or lower abdominal segments to decrease postoperative pain.
- Designed to prevent infections since separation of drainage systems, driver line systems or other devices is eliminated. The bands simply close above and below such items which exit at the slit(s) between them.
- Provides a steady pleasant support to a tender wound which reduces pain from bruising and swelling, improves blood circulation and oxygen supply to the surgical site to speed up healing of the wound.
- Provides a supportive pressure to the abdominal muscles to reduce pain during increased internal abdominal pressures when coughing, sneezing and straining.
- One or two of the bands can be opened for wound access while the other band(s) remain closed to maintain support and make wound care less painful.
- Provides support to the lumbar and abdominal regions which helps to get out of bed sooner, improves posture, breathing and mobility, and to walk around more easily and safer with the hands free.
- Converts easily from a 2 or 3 band support to a 1 or 2 band for additional comfort when sitting down.
- Convenient universal device holder available to arrange and attach drainage systems directly to the binder (sold separately).
- For more information on the product features, sizes and use, please refer to the product brochure or visit www.qualiteam.com

Sizing guide and order codes:

Supplied single packed with Instructions for Use

| | 1. Measure the fullest part of the abdomen, hips or waist. | Measurement 2: Less than 20 cm | Measurement 2: More than 20 cm |
|-------------------------------------|--|--------------------------------|--------------------------------|
| | 2. Measure the distance from breast bone to pubic bone. | Two band Height: 10 cm | Three band Height: 15 cm |
| Measurement 1 | Size | Order code | |
| 35 - 50 cm | Toddler | QLY2TOD | QLYTOD |
| 50 - 65 cm | Pediatric | QLY2PED | QLYPED |
| Custom sizes available upon request | | | |

| Accessories | |
|--|------------|
| Description | Order code |
| QualiPad Pediatric seatbelt pad | QTSEPAP |
| Universal device holder | DHG |
| Instructional plushy kit with QualiBelly Pediatric | LKY-QLY |

| Postoperative care kit | |
|--|------------------|
| Postoperative care kit (QualiBelly Advanced Ped. + QualiPad Ped. + Plushy kit supplied in washing net) | QTKIT + QLY size |



Free of charge:

Pediatric products are supplied with a plushy toy dog wearing a mini QualiBelly Advanced, along with color pencils and a coloring booklet illustrating "Instructions for Use" of QualiBelly Advanced.

Toddler products are supplied with the plushy toy dog only.

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.

For more information please contact Qualiteam at info@qualiteam.com
For orders email: orders@qualiteam.com

QUALIBelly


Post-partum, lumbar and abdominal binder



- Nurse designed compression binder ideal for postpartum and post-surgical wound support to prevent hernia and dehiscence after interventions such as Cesarean section, Liposuction, Abdominoplasty, and smaller interventions in the abdominal area.
- Provides a supportive pressure to the abdominal muscles which reduces pain during increased internal abdominal pressures when coughing, sneezing and straining.
- Provides a steady pleasant support to a tender wound which reduces pain from bruising and swelling, improves blood circulation and oxygen supply to the surgical site, which will speed up healing of the wound.
- Provides support to the lumbar and abdominal regions which helps to get out of bed sooner, improves posture, breathing and mobility, and to walk around more easily and safer with the hands free.
- Provides a comfortable compression to help loose skin shrink and reattach to the abdominal muscles, to get back into shape postpartum, post-cesarean, and after significant weight loss.
- For more information on the product features, sizes and use, please refer to the product brochure or visit www.qualiteam.com

Sizing guide and order codes:

Supplied single packed with Instructions for Use

|  <p>Measure the fullest part of the abdomen, hips or waist</p> | | |
|---|-------------|------------|
| Measurement | Size | Order code |
| 65 - 85 cm | Small | QLY20-S |
| 85 - 105 cm | Medium | QLY20-M |
| 105 - 125 cm | Large | QLY20-L |
| 125 - 145 cm | Extra Large | QLY20-XL |
| 145 - 165 cm | XX Large | QLY20-XXL |
| 165 - 185 cm | XXX Large | QLY20-XXXL |
| Custom sizes available upon request | | |

| Accessories | |
|-----------------------|------------|
| Description | Order code |
| QualiPad seatbelt pad | QTSEPA |

| Postoperative care kit | |
|---|-----------------------|
| Postoperative care kit (QualiBelly + QualiPad supplied in washing net) | QTKIT + QLY20 size |

Does not contain natural rubber latex (no NRL).

Washable at 30°C. Do not tumble dry.

Accessories:

CryoPouch

Non-drug adjunctive treatment tool to reduce pain.



- Two-pocket pouch with gel packs and fixation straps.
- Cool relief of incision pain in chest or abdomen, and of pain from bruised or broken ribs or sports injuries.
- Easy attachment to chest, abdominal binders, an arm or knee.
- Practical tool for healthcare staff to deliver non-drug pain relief as an adjunctive pain therapy.
- Hands-free cold treatment with increased convenience for patients and staff.
- For more information on the product features and use, please refer to the product brochure or visit www.callabyqualiteam.com

| Description | Order codes |
|---|-------------|
| CryoPouch two-pocket pouch with 2 cooling packs and fixation straps | CP1 |

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.

Accessories:

Reusable dual-purpose cool or heat packs

Gel packs for cool or heat treatment



- Reusable dual-purpose gel pack for cold treatment which relieves minor pain, prevents swelling and bruising, - or heat treatment which relieves minor muscle and joint pain, stiffness and muscle spasm.
- Special size made for insertion into CryoPouch.
- Also useful to reduce pain, bruising and inflammation after other injuries.
- Both cold and heat treatment should not be applied directly on the skin but always in a protective pouch
- Should not be used for more than 20 min periods. Can be repeated every 2 hours or as recommended by a health care professional.
- For more information on the product features and use, please refer to the product brochure or visit www.qualiteam.com

| Description | Order code |
|--------------------|------------|
| Gel pack 14 x 6 cm | CLD |

For more information please contact Qualiteam at info@qualiteam.com
For orders email: orders@qualiteam.com




Accessories:

AxillaPad

Perspiration pads for use with QualiBreath, QualiBra Advanced or Thor products



- Soft cotton pads to protect wound support products such as QualiBreath, QualiBra Advanced and Thor from excessive perspiration.
- Keeps the wound support products fresh for longer.
- Provides increased comfort under the armpits which might be needed when patients are bedridden for longer periods.
- Folds easily over the upper border of the wound supports and attaches with velcro.
- Axilla Pads are sold separately as well as included in postoperative care kits.
- For more information on the product features and use, please refer to the product brochure or visit www.qualiteam.com

| Version | Order codes |
|--|-------------|
|  AxillaPad perspiration pads | QBAX |
|  Pediatric AxillaPad perspiration pads | QBAXP |
|  Infant-Toddler AxillaPad perspiration pads | QBITAX |

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.



Accessories:

QUALIPad

Postoperative wound protection



- Nurse designed seat belt pad with a generous polyurethane filling. Protects surgical wounds, trauma areas in chest or abdominal region from the seat belt. Improves comfort for pregnant and obese persons, long distance and taxi drivers.
- Snaps easily on to the diagonal or lap section of the seat belt with Velcro straps and keeps the seat belt away from a sore area.
- Provides a buffer effect to protect a wound or trauma area at any impact pressure.
- Remains on the seat belt in the car.
- QualiPad is sold separately as well as included in postoperative care kits.
- For more information on the product features and use, please refer to the product brochure or visit www.qualiteam.com

| Version | Order codes |
|--|-------------|
|  QualiPad - seat belt pad | QTSEPA |
|  QualiPad Pediatric - seat belt pad | QTSEPAP |

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.

For more information please contact Qualiteam at info@qualiteam.com
For orders email: orders@qualiteam.com

Accessories:


Breast Cover

For females with breast size below cup D wearing QualiBreath.



- The breast cover is NOT A BRA but a cover for breast tissue appearing below the QualiBreath.
- Attaches with velcro at each upper corner to the inside of QualiBreath, and has a chest surrounding elastic band with a closure at the front left side.
- Adds comfort and privacy when women with smaller breast sizes use QualiBreath.
- For breast cups D and larger we recommend to use the QualiBra Advanced - breast and sternum support (see page 10-11)
- For more information on the product features and use, please refer to the product brochure or visit www.qualiteam.com

Supplied single packed with instructions for use

|  <p>Measure the circumference under the bust</p> | | |
|---|-------|-------------|
| Measurement | Size | Order codes |
| 75 - 110 cm | Small | COV-S |
| 110 - 145 cm | Large | COV-L |

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.


Accessories:

Universal device holder and line organiser

For organisation of drainage systems or other devices



- Drainage bottles or other devices can be secured and fastened by velcro in the holder.
- The holder is arranged to hang directly around one of the bands on QualiBelly Advanced, Calla Abdominal Compression Pro or the Thor.2 chest compression binder without any use of pins or tape.
- Drainage lines or other lines exiting from the abdomen or chest are easily secured by velcro on the holder.
- As the holder is attached directly on the compression binders, patients with drainage systems or other devices do not need to wait for staff to “liberate” them from bed. This helps to make patients more independent and to get active sooner.
- Available in green color for the QualiBelly Advanced abdominal binder, and in black color for the Calla abdominal compression Pro and the Thor.2 compression binder.
- For more information on the product features and use, please refer to the product brochure or visit www.qualiteam.com

| Versions | | Order codes |
|---|-------------------------------|-------------|
|  | Universal device holder green | DHG |

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.

For more information please contact Qualiteam at info@qualiteam.com
For orders email: orders@qualiteam.com

Accessories:

Netting Bag

for postoperative care kits



- All postoperative care kits are delivered in this convenient netting bag.
- The netting bag has a double function both as a container for the postoperative kit items and as a washing net to protect other garments when products with velcro closures are being washed.
- Space is provided on the netting bag to write the patient's name and department when hospital staff sends postoperative care products to a washing center.
- For more information on the product features and use, please refer to the product brochure or visit www.qualiteam.com

| Description | Order code |
|---|------------|
| Netting bag - for postoperative care kits | QT-NET |

QUALI CLIP

Disposable vascular bulldog clamp



- Atraumatic disposable bulldog clamp for temporary occlusion of smaller blood vessels
- Slightly serrated jaw pads for gentle, secure hold on the vessel
- Parallel clamping to optimize atraumatic occlusion
- Vivid color for easy distinguishing in the surgical field
- Radiopaque for effortless recognition
- Lightweight and Latex-free (no NRL)
- For more information on the product features and use, please refer to the product brochure or visit www.qualiteam.com

Sterile - individually packed

| Description | Units per box | Order code |
|------------------------------|---------------|------------|
| QualiClip disposable bulldog | 15 pcs | QT-011 |

Not available in the U.S.A.

QUALI *Maternity*




For advanced pregnancy care



- Nurse designed maternity support belt to relieve back and hip pain in the second and third trimester of pregnancy.
- Provides back support, improves posture and helps to prevent pelvic instability.
- Stabilizes the pelvis which helps to decrease Pelvic Girdle Pain (PGP) and to prevent separation of the pubic joint.
- Lifts and supports the expanding belly to protect the abdominal muscles and ligaments from over-stretching which may prevent or decrease stretch marks, varicoses and swelling of legs.
- Reduces pressure on the bladder to prevent potential incontinence.
- For more information on the product features, sizes and use, please refer to the product brochure or visit www.qualimaternity.com

Sizing guide and order codes:

Supplied single packed with Instructions for Use

| Measure around hips and below belly | | 10 cm band | 16 cm band |
|---|--------|---|---|
|  | |  |  |
| Measurement | Size | Order code | |
| 55-85 cm | Small | QM10S | QM16S |
| 85-115 cm | Medium | QM10M | QM16M |
| 115-145 cm | Large | QM10L | QM16L |
| Custom sizes available upon request | | | |

| Accessories | |
|-----------------------|------------|
| Description | Order code |
| QualiPad seatbelt pad | QTSEPA |

| Postoperative care kit | |
|---|--------------------|
| Postoperative care kit (QualiMaternity + QualiPad supplied in washing net) | QTKIT + QM size |

Does not contain natural rubber latex (no NRL).
Washable at 30°C. Do not tumble dry.

Gifts for patients:

T-shirts

Get well soon gift



- Flowers, plants, food, drinks and balloons are banned in many hospital wards as gifts to a patient after surgery.
- These awesome, high quality T-shirts of 100% cotton are great as a get-well-soon gift to cheer up a family member, friend or co-worker after surgery.
- Can be purchased separately or added to a postoperative care kit
- For more information on the product features and use, please visit www.qualiteam.com

| Lenght | Width | Size |
|--------|-------|-------------------|
| 74 cm | 51 cm | Medium |
| 76 cm | 56 cm | Large |
| 79 cm | 61 cm | Extra Large |
| 80 cm | 67 cm | Extra Extra Large |
| 84 cm | 71 cm | 3 x Extra Large |
| 86 cm | 76 cm | 4 x Extra Large |
| 89 cm | 81 cm | 5 x Extra Large |

Can only be ordered on www.qualiteam.com where designs and size can be selected.


Gifts for patients:

Silk scarves

Get well soon gift




- Get-well-soon gift
- Flowers, plants, food, drinks and balloons are banned in many hospital wards as gifts to a patient after surgery.
- This beautiful Italian made, 100% silk scarf is a wonderful get-well-soon gift to cheer up a family member, friend or co-worker, and to celebrate a successful surgery.
- Can be purchased separately or added to a postoperative care kit
- For more information on the product features and use, please visit www.qualiteam.com

|  | Description | Order codes |
|---|---------------------------------|-------------|
| | Qualiteam silk scarf 70 x 70 cm | SCA70 |
| | Qualiteam silk scarf 90 x 90 cm | SCA90 |

Instructional plushy toy



- [illegible]

| Versions | Order codes |
|---|--|
|  <p data-bbox="388 1016 504 1037">Lucky - Instructional plushy toy</p> | <p data-bbox="555 1016 625 1037">LKY_01</p> |
|  <p data-bbox="388 1143 504 1174">Lucky - Instructional plushy toy wearing QualiBreatH</p> | <p data-bbox="555 1143 625 1164">LKY-QB</p> |
|  <p data-bbox="388 1269 504 1300">Lucky - Instructional plushy toy wearing QualiBelly Advanced</p> | <p data-bbox="555 1269 625 1291">LKY-QLY</p> |



Copyright © Qualiteam Group Ltd 2021
All rights reserved