

D'VILL™

Long-length, kink-resistant hemostasis
introducer for introduction of large
diameter interventional devices

Hemostasis valve
designed to minimize
blood loss

**Side port and
3-way stopcock**
for convenient flushing


Braided sheath
for optimal flexibility
and kink resistance

**Radiopaque
marker band**
for easy imaging

**Available in
30 cm, 65 cm,
and 85 cm
lengths**

FOR YOUR MOST COMPLEX CASES

The D'VILL™ is made of an innovative braided material for enhanced flexibility and kink resistance. Available in 12F and 14F sizes, the D'VILL™ is compatible with CP Stent® and BIB® catheters.



Indications for use:

Recommended for introduction of balloons, catheters, and other diagnostic and interventional devices. Refer to the D'VILL™ Instructions for Use for relevant warnings, precautions, complications, and contraindications.

D'VILL™ Introducer

Long-length, kink-resistant
hemostasis introducer



For Introduction of Large Diameter Interventional Devices

Compatible with CP Stent® and BIB® catheters

D'VILL™ Specifications			
Catalog No.	Introducer Size (FR)	Usable Length (CM)	Guidewire (inches)
DV1230	12	30	0.035
DV1430	14	30	0.035
DV1265	12	65	0.035
DV1465	14	65	0.035
DV1285	12	85	0.035
DV1485	14	85	0.035

- Each D'VILL™ set includes (1) one sheath and (1) one dilator

The D'VILL™ is available worldwide. Please contact us or your local distributor to place an order.

This device is subject to individual country regulations regarding the importation and sale of this product



G. van Wageningen B.V. European Representative

info@heartmedical.nl
Hallenweg 40
5683 CT Best
The Netherlands
Tel: +31-499-377388
Fax: +31-499-377456

NuMED | 2880 Main Street | Hopkinton, NY 12965 USA
Tel: 315.328.4491 | Fax: 315.328.4941 | www.numedforchildren.com

