



The **Hemostasis Introducer** is recommended for the use with the pfm medical family of Nit-Occlud® occlusion systems.

Details

- ▶ **Hemostatic Valve** for all french sizes up to 18F for easy loading and prevention of backbleeding
- ▶ **Mesh braiding** for excellent kink resistance and pushability
- ▶ **PTFE inner layer** enables smooth and easy delivery of occluders and other devices
- ▶ **Smooth transition** between dilatator and tip for smooth advancement
- ▶ **Platinum iridium Marker band** for enhanced visibility under fluoroscopy

Ordering Information

REF	Sheath Diameter	Curve	Sheath Length	Connection	Braiding	PU
95040851*	4F	45°	85 cm	Hemostatic valve	Mesh	1
95050851*	5F	45°	85 cm	Hemostatic valve	Mesh	1
95060851	6F	45°	85 cm	Hemostatic valve	Mesh	1
95070851*	7F	45°	85 cm	Hemostatic valve	Mesh	1
95080851	8F	45°	85 cm	Hemostatic valve	Mesh	1
95090851	9F	45°	85 cm	Hemostatic valve	Mesh	1
95100851	10F	45°	85 cm	Hemostatic valve	Mesh	1
95120851	12F	45°	85 cm	Hemostatic valve	Mesh	1
95140851	14F	45°	85 cm	Hemostatic valve	Mesh	1
95160851*	16F	45°	85 cm	Hemostatic valve	Mesh	1
95180851*	18F	45°	85 cm	Hemostatic valve	Mesh	1

* Only available upon request

Contact

Should you have any questions our Customer Solutions Team will be glad to assist you.

✉ service@pfmmedical.com

☎ +49 (0)2236 9641-221

☎ +49 (0)2236 9641-20

Distributor
pfm medical ag
Wankelstraße 60
50996 Köln, Germany

Certified according to
DIN EN ISO 13485

Manufacturer
APT Medical Inc.
No.009, Xiangxiang Economic, Development Zone,
Xiangxiang Hunan 411400 China

CE 0120